Nhe I



5'G▼CTAGC····3' 3'CGATC▲G···5'

Nhel is a restriction enzyme purified from *Neisseria mucosa heildelbergensis* (ATCC 25999).

Catalogue No 146-1, 1000 U

146-2, 3x1000 U

Concentration 10-12u/μl and 40-

60u/µl*

Reagents supplied: 10x A and 10x K buffer

Unit substrate: Lambda DNA (*Hind*III digest).

Unit calculation assay conditions: 50 mM potassium acetate, 20 mM Tris-acetate (pH 7.9 @ 25°C), 10 mM magnesium acetate, 1 mM dithiothreitol, 100 μ g/ml BSA. Incubate at 37°C.

Absence of contaminants: 80 units of Nhel do not produce any unspecific cleavage products after 16 hrs incubation with 1 μ g of λ DNA/Hind III digest at 37°C. After 100-fold overdigestion with Nhel, greater than 98% of the DNA fragments can be ligated and recut.

Storage buffer: 200 mM NaCl, 10 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM dithiothreitol, 0.15% Triton X-100, 200 μ g/ml BSA and 50% glycerol. Store at -20°C.

Heat inactivation: 65°C for 20 minutes...

Methylation Sensitivity:

dam methylation: Not sensitive dcm methylation: Not sensitive

CpG methylation: Blocked by some combinations of overlapping

Note: Activity inhibited by salt concentrations >100mM. Cleaves to leave a 5' CTAG extension which can be efficiently ligated to DNA fragments generated by *AvrII*, *SpeI*, or *XbaI*.

Star activity: Low salt, high glycerol (>5%) concentrations, pH >8.0 or large excess of the enzyme may result in star activity.

Percent Activity in MINOTECH Buffers

L	М	Н	SH	Α	К
100	50-75	0-20	<10	100	100

General reaction mixture:

10U Nhel	1μΙ
10x A or K buffer *	2μΙ
DNA substrate	<1µg
Sterile ultrapure water	Up to 20 μl
Incubate for 15 r	nin at 37°C

^{*}In the case of A buffer we recommend the addition of BSA to a final concentration of $100 \mu g/ml$.

Frequency of Cutting

λ	Ad-2	Фх174	pUC18	M13mp18	pBR322
1	4	0	0	0	1



Lambda DNA 0.7 % agarose



^{*}Add an H to cat.# to order the high concentration